## FORCED AIR COOLING FIELD TRIAL RESULTS

## Watermelon



A precooling experiment was conducted to test the performance of a small scale (countertop sized) forced air cooler (FAC) in parallel with product cooled by room cooling. Freshly harvested watermelon (123 lbs)

was placed in a CoolBot controlled walk-in cooler set at 46 °F. In the case of room cooling the product temperature is reduced as a result of simply being in the room with cool air circulating around it. The same mass of product was placed in a small forced air cooler (FAC) that included a base, frame, suction fan, plenum, and plastic tarp with one end open to direct the cold room air over the product packed inside the crates. The ambient temperature of the cooler and the pulp temperature of the produce cooled using each method was monitored over time to determine and compare the precooling rate. The product started at 96 °F and, over the course of an hour, dropped 3 °F by room cooling and 7 °F by forced air cooling.





| FAC / RC Ratio ("FAC is times faster") |             | 1.6               |
|--|-------------|-------------------|
| Time to 7/8 Temp (Hours), estimated    | 29          | 14                |
| Observed Cooling Rate (degF/min)       | 0.10        | 0.17              |
| Temp @ Test End [2hr 35min], (°F)      | 79          | 72                |
| Temp @ 60min (°F)                      | 92          | 91                |
| Temp @ 20min (°F)                      | 94          | 97                |
| Starting Temp (°F)                     | 95          | 98                |
| Melons (123.4lbs)                      | Room Cooled | Forced Air Cooled |
| Cooler Set point (°F) 46               |             |                   |
|  |             |                   |

A standard measure of precooling rate is the time required to bring the product down % of the way to the target storage temperature. This is called "% time". Based on this test, when starting at 96 °F, it was estimated that the % time for forced air cooling was 14 hours and for room cooling it was 29 hours. These results show that it takes 2.1 times longer to room cool watermelon when compared to FAC (or FAC is 1.6 times faster).

## **Acknowledgments**

Funding for this publication was made possible, in part, by the USDA NE SARE program under grant #LNE16-347. Thanks to Clear Brook Farm for participation in this trial.





This publication is available at go.uvm.edu/factrial

Andrew S. Chamberlin and Christopher W. Callahan



